

PERFORMANCE ANALYSIS OF HIGH CAPACITY MASS COMPARATOR- AX64004

Xiaoping Ren/Presenter¹, Jian Wang¹, Ruilin Zhong¹, Changqing Cai¹, and Qingxiong Ren²

(1. Division of Mechanics and Acoustic, National Institute of Metrology, Beijing 100013, China;
(2. Shanxi Province Institute of Metrology, Taiyuan 030002, China)

Abstract: The calibration of heavy mass standards is a difficult, demanding, and time consuming process. Loading and unloading of heavy mass standards is strenuous work as well as placing of the weights in the center of the weighing position by human. A high capacity automatic comparator with readability of 0.1 mg has become available in National Institute of Metrology (NIM). Performance analysis of this mass comparator is shown in the paper.

Keywords: Mass Dissemination; Automatic Mass Comparator; Performance analysis.

1. INTRODUCTION

Recently, a new generation of computer controlled comparators-AX64004, particularly designed for the calibration of high capacity weight, has become available. This equipment can automatically executes all steps involved in disseminations, including the moving of the weights and the necessary comparisons, requiring neither manual intervention by the operator, nor custom auxiliary weights [1].

The requirements of high-capacity mass dissemination to calibrate weights and weight sets in the nominal range of a few kilograms up to 60 kg are first discussed by Reichmuth and Richard in [2]. Now there are several national metrology institutes and independent commercial mass metrology laboratories equipping this comparator. For example: troemner is an independent calibration laboratory in the United States that can provide calibrations with this type of balance [3]. Primary Standards Laboratory of Sandia National Laboratories also maintains this kind of balance[4]. In Klodian Dhoska's master thesis [5], he used this kind of comparator to get experience in the field of high accuracy mass measurements. In [6], the calibration of multiples of 1 kg is carried out using a fully automatic mass comparator with maximum capacity of 64 kg and resolution of 0.1 mg. In Jan.26th2010, China's first AX64004 was settled in Shanghai verification and testing technology institute [7].

Form all these documents and websites ,we can see that the requirements of high-capacity mass dissemination are very popular, like Estonia[8], 50 kg comparison of Euromet from 2004 to 2007 [9] and CCM.M-K3.1 comparison of 50 kg mass in 2009 to 2010[10].

2. AUTOMATIC MASS COMPARATOR-AX64004

AX64004 is an automatic mass comparator which is based on computer operation, has four positions for weight alternator [6]. This electronic mass comparator is from Mettler Toledo, Switzerland, and specially used for large precision weights from 10 kg to 60 kg (weighing capacity is 64 kg, with remarkable high resolution of 0.1 mg and a typical standard deviation of less than 0.4 mg [5].

The AX64004 comparator weighing system comprises [5]: weight handler; balance controller; handler controller; temperature sensors; 4 glass cylinders (individual draft shields, designed for conventional 50 kg normal significantly with larger dimensions); cables connecting to weighing unit; computer control software; crane used for lifting 50 kg weights. That is all for fully automatic performance and reporting of comparative weightings. Besides, it has a modern and fast weighing system with switching weights (250g-steps). Comparator's substitution weight set is shown at the bottom of equipment. The weight partition provides combinations of substitution loads from 0 to 64 kg in increments of 250 g. With the weighing cell's electrical range of 260 g, this allows for an overlapping weighing range from 0 to 64 kg. Also kept with this set is a calibration weight of 250 g [2].



Figure 1 AX64004 comparator installed in NIM and crane used for lifting weight

It features a four-place turntable for automatic weight ex-change, which is shown in the foreground in Fig1. The weighing platforms, equipped with individual draft covers, offer enough space to accommodate weight combinations according to any weighing scheme. METAS has given out an example on large object calibration without draft shield

[2]. Besides, the weights self-centre mechanism can reduce the eccentric load errors.

The comparator's display unit and a laptop PC to control the comparator, is visible on the right in the background. Computer control software with NIM's logo shown in Fig 2 is designed by National Institute of Metrology, China (NIM) and Mettler Toledo Corporation (MT).



Figure 2. Computer control software designed by NIM and MT

3. PERFORMANCE ANALYSIS OF AX64004 MASS COMPARATOR

On Apr.29th 2013(Apr 29th 09:22:37 to May 3rd 13:14:34, and it is holiday from Apr 29th to May 1st), an equipment test was generated by AX64004 according to running an all comparison weighing processes consisting of 6 series of 10 A-B-A cycles weighing. And the running data was analysed. The comparisons are P1-P1, P2-P1, P3-P1, P4-P1, P2-P2, P3-P2, P4-P2, P3-P3, P4-P3, P4-P4, because the weighing scheme is A-B-A cycles, so we make all the comparison cycle to the same start point. Weighing data of P1, P2, P3 and P4 position are shown in figure 3. These measurement record obtained from 6 series of 10 comparisons between the same positions used one 50 kg mass standards.







Figure 3. Measurement data from the same weighing positions

Table 1 is the mean (conventional) standard deviation of the same position, calculated from groups of 10 A-B-A comparisons (mark by circle, multiply sign, triangle, plus, square, star). Here in the weighing cycles, "A" represents weighing the reference weight and "B" represents weighing the test weight. The "A-B-A" cycle is normally used when calibrating class E and F weights automatically.

The maximum value of each position is 0.11, 0.10, 0.09 and 0.21. All the values are better than the parameter shown in [6][2](0.4/0.2 mg), [5](0.29 mg).

Table 1 Standard Deviation of the same weighing position

Series No.	Standard Deviation (mg)			
	P1	P2	P3	P4
1	0.11	0.10	0.04	0.04
2	0.04	0.04	0.05	0.06
3	0.05	0.02	0.05	0.06
4	0.04	0.03	0.07	0.03
5	0.06	0.09	0.07	0.09
6	0.05	0.04	0.09	0.21




Figure4. Measurement standard deviation obtained from 60 series comparisons between two weight positions

Figure 4 shows the measurement standard deviation obtained from 60 series comparisons between two weight positions. We can see that during P3 and P2 comparison, the standard deviations of first and second series are 0.50 and 0.52 respectively, which are more higher than 0.4 mg. Maybe it is due to a transient at the beginning of the measurement, probably caused by acclimatization. Another reason for the data is eccentric load error. Although it has self-centre mechanism, the 50 kg weight is too heavy and it is hard to centre the weight at only a few centring process.

There is another high standard deviation (0.45 mg) happens to Series 6 between P4 and P3 position. It was almost the last comparison of the whole test, and holiday was finished. Human and environment reason can result to this value. Besides, the other measurement standard deviation almost remains smaller than 0.25 mg for the remaining measurements.

4. CONCLUSIONS

The repeatability of AX64004 comparator is always less than 0.4 mg (some times even better than 0.2 mg), which demonstrates that required environmental conditions indispensable for the running performance of this instrument.

Only after satisfactory repeatability results are given, other characteristics can be tested and usefully determined. Besides, special care should to be given to suppress

eccentric load error. For example, 50 kg kilogram weight is so heavy and big that, it cannot always be easily centred.

5. ACKNOWLEDGEMENT

The author wish to thank for the grant of National Science and Technology Support Program (2011BAK15B06); Grant of Special-Funded Program on National Key Scientific Instruments and Equipment Development (2012YQ090208), and China Postdoctoral Science Foundation (2013M541003).

6. REFERENCES

- [1] Arthur Reichmuth. A New Mass Comparator Generation for the Automatic Calibration of Weight Sets. Proceedings of the 17th International Conference on Force, Mass, Torque and Pressure Measurements, IMEKO TC3, 17-21 Sept. 2001, 310-319.
- [2] Arthur Reichmuth, Philippe Richard. High-capacity mass dissemination with four place mass comparator. San Diego 2002 Workshop and Symposium, <http://www.ncsli.org/c/a/h/nl/2003-01-NLa.pdf>
- [3] Troemner Calibration Services–Mass. <http://www.troemner.com/pdf/deskref/TroemnerCalibrationServicesMass.pdf>
- [4] Sandia National Laboratories. Primary Standards Laboratory Metrology Program. http://www.sandia.gov/psl/MSP_Fact_Sheets_082708/MASS_FORCE_factSheet2008.pdf
- [5] Klodian Dhoska. Systematic effects in automated mass measurements. Master thesis of University of Tartu, 2010.
- [6] Goutam Mandal, C. K. Gopan. Improvement of measurement capabilities of multiple mass standards of 1 kg at NPLI. IMEKO 2010 TC3 Conference Metrology in Modern Context 2010, Pattaya, Chonburi, Thailand, 2010, 87-90.
- [7] http://news.frbiz.com/china_39:s_first_fully_automatic-180334.html
- [8] Lars Nielsen. Report of EUROMET Mass & Related Quantities TC Meeting, 2005, 117. http://resource.npl.co.uk/docs/science_technology/mass_force_pressure/euromet/2005.pdf
- [9] André Evenstad JV, Chris Mitsas EIM. Report of Euromet 83 2: 50 kg comparison. http://www.mikes.fi/documents/upload/j5_2007.pdf
- [10] G Mann, T Madec, P A Meury. Report of CCM.M-K3.1 comparison/50 kg mass. http://kcdb.bipm.org/appendixb/appbreports/ccm.m-k3/ccm.m-k3.1_final_report.pdf